



TEST KIT INSTRUCTIONS

Half fill the plastic test tube with water to be tested. Add 1 drop of solution from dropper bottle to the test tube, fit test tube cap and shake solution vigorously.

A strong lather with a clear solution indicates water is fully 'soft'. This is calibrated to occur at a hardness of less than 14.3mg/l CaCO_3 (1° clarke).

No lather and a milky solution indicates water is 'hard' and softener needs regeneration or attention. Continuing to add drops until a lather is achieved will give guidance on severity of water hardness; i.e. 1 drop = fully 'soft', 8 to 12 drops = fully 'hard'.

Hard water is typically 286mg/l CaCO_3 (20° clarke).

www.TheWaterSoftenerCentre.com



Tablet Salt & Block Salt
collection discounts and delivery service

FOR YOUR NEXT SALT :-

Call
Click
Collect

Tel 01672 518040
www.TheWaterSoftenerCentre.com
Ensign, Unit 6 Hertford Court
Marlborough Business Park
Marlborough SN8 4AW

Water Softener Sales,
Rentals & Installations

Water Softener Service,
Repairs & Trade-in's

Mon - Fri 8:00am - 5:00pm
Sat 9:00am - 1:00pm
Sun closed

ENSIGN 

The Water Softener Centre

www.TheWaterSoftenerCentre.com

Water Softener Specialists Since 1982

Tel 01672 518040